

## Experiments to Learn Primer Sequence Bias of DNA Polymerase

### Overview

The purpose of this set of experiments is to learn the DNA polymerase (DNAPol) priming efficiency of all 7mers. This will be done with primers that contain degenerate nucleotides (e.g. Random\_P5 = GATCTCCGAGTTGCNNNNNNN), and the goal is to evaluate the priming efficiency of all  $4^7$  combinations of NNNNNNN. Unfortunately, strong primer hairpins will form for certain NNNNNNN, and these NNNNNNN will not prime DNAPol and will falsely indicate very poor priming efficiency. To overcome this problem, there are two sets of experiments ("DNAPol Experiments 1" and "DNAPol Experiments 2", below) distinguished by having different "random primers" that have been designed such that the hairpins that form for one random primer will not form in the other random primer. The combined results of the two experiments will be used to quantify the DNAPol priming efficiency of all 7mers.

### Primer Sequences

bkgrnd\_P5 = CGACGCTCTTCCGATCTGGTTTGTCTACGAGGAGG

bkgrnd\_P7 = GTGTGCTCTTCCGATCCCCCGTCAAGTAGGAGG

Common\_P5 = CGACGCTCTTCCGATCTACTCTGTCCGATGTCCTGCCTTGCTTTTGGTCCG

Common\_P7 = GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTTGGTCCG

Random\_P5 = GATCTCCGAGTTGCNNNNNNN

Random\_P7 = CGATCCGCTCAACTNNNNNNN

Univ\_P5+ = CGACGCTCTTCCGATCTCCGAG

Univ\_P7+ = GTGTGCTCTTCCGATCCGCTC

### Primer Tm's

Random\_P5 = 60-70C (for NNNNNNN = AAAAAAA,...,GGGGGGG)

Random\_P7 = 60-70C (for NNNNNNN = AAAAAAA,...,GGGGGGG)

bkgrnd\_P5 = 58C (for the initial hybridization region)

bkgrnd\_P7 = 60C (for the initial hybridization region)

Univ\_P5+ = 35C (for partial hybridization)

= 67C (for full hybridization)

Univ\_P7+ = 42C (for partial hybridization)

= 67C (for full hybridization)

Universal primers (Used to construct Univ\_P5+ and Univ\_P7+; included here to confirm correct sequences use)

>UDTD5

ACACTCTTTCCCTACACGACGCTCTTCCGATCT

P5\_1                   CGACGCTCTTCCGAT

>UDTD7

GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC

P7\_1                   GTGTGCTCTTCCGAT

Protocol Outline for all experiments

- 1) For each of the experiments below, do only 1 cycle of PCR with specified primers.
- 2) PCR cleanup with Ampure beads
- 3) 15 cycles of PCR with the augmented universal primers Univ\_P5+ and Univ\_P7+.
- 4) PCR with Illumina sequencing adapters.

## DNAPol Experiment 1

### Experiments to Perform

#### 1.1) DNAPol Experiment 1.Back

- 1 cycle PCR with DNAPol\_template\_1 as input DNA and the two primers bkgrnd\_P5 & Common\_P7
- PCR Cleanup
- Standard PCR with the primers Univ\_P5+ and Univ\_P7+

#### 1.2) DNAPol Experiment 1.Bias

- 1 cycle PCR with DNAPol\_template\_1 as input DNA and the two primers bkgrnd\_P5 & Common\_P7
- PCR cleanup
- Standard PCR with the primers Univ\_P5+ and Univ\_P7+

### Template and Primers

Common_P7		Random_P5	
5'-GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCG		NNNNNNCGTTGAGCCTCTAG-5'	
	5'-CCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNGCAACTCGGAGATCCCTCCTCGTAGACAAACC-3'		bkgrnd_P5
DNAPol_template_1			GGAGGAGCATCTGTTTGGTCTAGCCTTCTCGCAGC-5'
			GAGCCTCTAGCCTTCTCGCAGC-5'
			Univ_P5+

### DNAPol Experiment 1.Back

"Background" product after PCR with bkgrnd\_P5 & Common\_P7 primers is 135bp:

5'-GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNGCAACTCGGAGATCCCTCCTCGTAGACAAACCAGATCGGAAGAGCGTCG-3'  
3'-CACACGAGAAGGCTAGTGAGACAGGCTACAGGACGGAACGAAAACCAGGCCAGANNNNNNNACATTAGGACCTCNNNNNNNNNCGTTGAGCCTCTAGGGAGGAGCATCTGTTTGGTCTAGCCTTCTCGCAGC-5'

### DNAPol Experiment 1.Bias

"Background" product after PCR with bkgrnd\_P5 & Common\_P7 primers is 100bp:

5'-GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNGCAACTCGGAGATC-3'  
3'-CACACGAGAAGGCTAGTGAGACAGGCTACAGGACGGAACGAAAACCAGGCCAGANNNNNNNACATTAGGACCTCNNNNNNNNNCGTTGAGCCTCTAG-5'

## DNAPol Experiment 2

### Experiments to Perform

#### 2.1) DNAPol Experiment 2.Back

- PCR with DNAPol\_template\_2 as input DNA and the two primers bkgrnd\_P7 & Common\_P5
- PCR cleanup
- Standard PCR with the primers Univ\_P5+ and Univ\_P7+

#### 2.2) DNAPol Experiment 2.Bias

- PCR with DNAPols\_template\_2 as input DNA and the two primers Random\_P7 & Common\_P5
- PCR cleanup
- Standard PCR with the primers Univ\_P5+ and Univ\_P7+

### Template and Primers

Common_P5		Random_P7	
5'-CGACGCTCTCCGATCTACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCG		NNNNNNNTCAACTCGCCTAGC-5'	
	5'-CCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNNNAGTTGAGCGGATCGCCTCCTACTTGACGGGG-3'		bkgrnd_P7
	DNAPol_template_2		GGAGGATGAAGTCCCCCTAGCCTTCTCGTGTG-5'
		CTCGCCTAGCCTTCTCGTGTG-5'	
		Univ_P7+	

### DNAPol Experiment 2.Back

"Background" product after PCR with bkgrnd\_P7 & Common\_P5 primers is 137bp:

5'-CGACGCTCTCCGATCTACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNNNAGTTGAGCGGATCGCCTCCTACTTGACGGGGGATCGGAAGAGCACAC-3'  
3'-GCTGCGAGAAGGCTAGATGAGACAGGCTACAGGACGGAACGAAAACCAGGCCAGANNNNNNNNACATTAGGACCTCNNNNNNNNNNNNNTCAACTCGCCTAGCGGAGGATGAAGTCCCCCTAGCCTTCTCGTGTG-5'

### DNAPol Experiment 2.Bias

"Bias" product after 1 cycle of PCR with Random\_P7 & Common\_P5 primers is 104bp:

5'-CGACGCTCTCCGATCTACTCTGTCCGATGTCCTGCCTTGCTTTGGTCCGGGTCTNNNNNNNTGTAATCCTGGAGNNNNNNNNNNNAGTTGAGCGGATCG-3'  
3'-GCTGCGAGAAGGCTAGATGAGACAGGCTACAGGACGGAACGAAAACCAGGCCAGANNNNNNNNACATTAGGACCTCNNNNNNNNNNNNNTCAACTCGCCTAGC-5'