

## Experiments to Learn Primer Sequence Bias of Reverse Transcriptase

### Overview

The purpose of this set of experiments is to learn the reverse transcriptase priming efficiency of all 7mers. This will be done with primers that contain degenerate nucleotides (e.g. Random\_P5 = GATCTCCGAGTTGCNNNNNNN), and the goal is to evaluate the priming efficiency of all  $4^7$  combinations of NNNNNNN. Unfortunately, strong primer hairpins will form for certain NNNNNNN, and these NNNNNNN will not prime RT and will falsely indicate very poor RT priming efficiency. To overcome this problem, there are two sets of experiments ("RT Experiments 1" and "RT Experiments 2", below) distinguished by having different "random primers" that have been designed such that the hairpins that form for one random primer will not form in the other random primer. The combined results of the two experiments will be used to quantify the RT priming efficiency of all 7mers.

### Primer Sequences

Random\_P5 = GATCTCCGAGTTGCNNNNNNN

Random\_P7 = CGATCCGCTCAACTNNNNNNN

Common\_P5 = CGACGCTCTTCCGATCT ACTCTGTCCGATGT CCTGCCTTGCTTTTGGTCCG

Common\_P7 = GTGTGCTCTTCCGATC ACTCTGTCCGATGT CCTGCCTTGCTTTTGGTCCG

Univ\_P5+ = CGACGCTCTTCCGATCTCCGAG

Univ\_P7+ = GTGTGCTCTTCCGATCCGCTC

### Primer Tm's

Random\_P5 = 60-70C (for NNNNNNN = AAAAAA,...,GGGGGGG)

Random\_P7 = 60-70C (for NNNNNNN = AAAAAA,...,GGGGGGG)

Common\_P5 = 60C

Common\_P7 = 60C

Univ\_P5+ = 30C (for partial hybridization)

62C (for complete hybridization)

Univ\_P7+ = 38C (for partial hybridization)

= 60.5C (for complete hybridization)

Universal primers (Used to construct Univ\_P5+ and Univ\_P7+; included here to confirm correct sequences used)

>UDTD5

ACACTCTTTCCCTACACGAGCTCTTCCGATCT

Universal\_P5 CGACGCTCTTCCGAT

>UDTD7

GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC

Universal\_P7 GTGTGCTCTTCCGAT

## RT Experiments 1

## Experiments to Perform

### 1.1) RT\_Experiment\_1.Back:

Standard PCR using RT\_template\_1 as input DNA and the two primers Univ\_P5+ and Common\_P7

### 1.2) RT\_Experiment\_1.Bias

RT-qSeq using RT\_template\_1 as "input RNA" and the two primers Random\_P5 and Common\_P7

## Template and Primers

Common\_P7 Random\_P5  
5'-GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTTGGTCCG NNNNNNNCGTGTGAGCCTCTAG-5'  
5'-CCTGCCTTGCTTTTGGTCCGGGTCTNNNNNTGTAATCCTGGAG NNNNNNNNNNNNN GCAACTCGGAGATC-3'  
RT\_template\_1 GAGCCTCTAGCCTTCTCTGCAGC-5'  
Univ\_P5+

Note: RNA bases (ribonucleotides) are indicated in red.

POTENTIAL PROBLEM: RNase treatment in the RT-qSeq for RT\_Experiment\_1.Bias would only partially degrade the "input RNA" RT\_template\_1. This could cause problems. Possible solutions: 1) include other RNA bases so that template is degraded to smaller pieces that would be removed in PCR cleanup; 2) no RNase treatment.

RT Experiment 1.Back

"Background" product after PCR with Univ P5+ & Common P7 primers is 112bp:

5' -GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGCTTTTGGTCCGGGTCTNNNNNTGTAATCCTGGAGNNNNNNNNNNNGCAACTCGGAGATCGGAAGAGCGTCG-3'  
3' -CACACGAGAAGGCTAGTGAGACAGGCTACAGGACGGAACGAAAACAGGCCAGANNNNNACATTAGGACCTCNNNNNNNNNNNNCGTTGAGCCTCTAGCCTTCTCGCAGC-5'

## RT Experiment 1.Bias

"Bias" product after "Qiagen 1-Step" step of RT-qSeq with Random P5 & Common P7 primers is 100bp:

5'-GTGTGCTCTTCCGATCACTCTGTCCGATGTCCTGCCTTGTCTTTTGGTCCGGGTCTNNNNNTGTAATCCTGGAGNNNNNNNNNNNGCAACTCGGAGATC-3'  
3'-CACACGAGAAGGCTAGTGAGACAGGCTACAGGACGGAACGAAAAACCAGGCCAGANNNNNACATTAGGACCTCNNNNNNNNNNNNCGTTGAGCCTCTAG-5'

## Experiments to Perform

Standard PCR using RT template 2 as input DNA and the two primers Univ P7+ and Common P5

RT-qSeq using RT template 2 as "input RNA" and the two primers Random P7 and Common P5

Common\_P5 Random\_P7

5'-CGACGCTCTTCCGATCTACTCTGTCCGATGTCCTGCCTTGCCTTTTGGTCCG NNNNNNNNTCAACTCGCCTAGC-5'

5'-CCTGCCTTGCCTTTTGGTCCGGGTCTNNNNNTGTAATCCTGGAG NNNNNNNNNNNNAGTTGAGCGGATCG-3'

RT\_template\_2 Univ\_P7+

CTCGCCTAGCCTTCTCGTGTG-5'

POTENTIAL PROBLEM: RNase treatment in the RT-qSeq for RT\_Experiment\_1.Bias would only partially degrade the "input RNA" RT\_template\_1. This could cause problems. Possible solutions: 1) include other RNA bases so that template is degraded to smaller pieces that would be removed in PCR cleanup; 2) no RNase treatment.

"Background" product after PCR with Univ P7+ & Common P5 primers is 112bp:

5'-CGACGCTCTTCGATCTACTCTGTCCGATGTCTGCCTTGCTTTTGGTCCGGGTCTNNNNNTGTAATCCTGGAGNNNNNNNNNNAGTTGAGCGGATCGGAAGAGCACAC-3'  
3'-GCTGCGAGAAGGCTAGATGAGACAGGCTACAGGACGGAACGAAACCAGGCCAGANNNNNACATTAGGACCTCNNNNNNNNNNNNNTCAACTCGCCTAGCCTTCTCGTGTG-5'

"Bias" product after "Qiagen 1-Step" step of RT-qSeq with Random P7 & Common P5 primers is 101bp:

5' -CGAGCGCTCTTCGATCTACTCTGTCCGATGTCCTGCCTTGTCTTTGGTCCGGGCTNNNNNTGTAATCTGGAGNNNNNNNNNNNAGTTGAGCGGATCG-3'  
3' -GCTGCGAGAAGGCTAGATGAGACAGGCTACAGGACGGAACGAAAACAGGCCAGANNNNNACATTAGGACCTCNNNNNNNNNNNTCAACTCGCCTAGC-5'